

KLINIPRESS - Nonwoven Compress Sterile

Description

Nonwoven swab for cleansing and covering light to non-exuding wounds. Soft and comfortable. Medical device. Sterile if pouch and seal are intact.

Product information

In compliance with: MDD (EU) 93/42/EEC Class Is



Article overview

Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
3017044	KLINION KLINIPRESS NONWOVEN COMPRESS STERILE 5 X 5 CM 4 PLY PACKED PER 100 X 1 PCS REF 175001 *S*	100 pcs.	2000 pcs.			175001
3017045	KLINION KLINIPRESS NONWOVEN COMPRESS STERILE 5 X 5 CM 4 PLY PACKED PER 50 X 2 PCS REF 175002 *S*	100 pcs.	2000 pcs.			175002
3017046	KLINION KLINIPRESS NONWOVEN COMPRESS STERILE 5 X 5 CM 4 PLY PACKED PER 20 X 5 PCS REF 175003 *S*	100 pcs.	3000 pcs.			175003
3017070	KLINION GAAS NONWOVEN KOMPRES EXTRA 5X5CM 8L *S* STERIEL 175101	100 pcs.	2000 pcs.			175101
3016479	KLINION KLINIPRESS NON WOVEN COMPRESS STERILE 5 X 5 CM PACKED PER 5 PCS 4 PLY REF 40512070 *S*	175 pcs.	5250 pcs.			40512070
3017048	KLINION KLINIPRESS NONWOVEN COMPRESS STERILE 7.50 X 7.50 CM 4 PLY PACKED PER 50 X 2 PCS REF 175008 *S*	100 pcs.	2000 pcs.			175008
3017049	KLINION KLINIPRESS NONWOVEN COMPRESS STERILE 7.50 X 7.50 CM 4 PLY PACKED PER 20 X 5 PCS REF 175009 *S*	100 pcs.	2000 pcs.			175009
3016480	KLINION KLINIPRESS NON WOVEN COMPRESS STERILE 7.50 X 7.50 CM PACKED PER 5 PCS 4 PLY REF 40512071 *S*	200 pcs.	3000 pcs.			40512071
3017052	KLINION KLINIPRESS NONWOVEN COMPRESS STERILE CONSUMER PACKAGE 10 X 10 CM 4 PLY PACKED PER 10 X 1 PCS REF 175029 *S*	10 pcs.	200 pcs.			175029

Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
3017053	KLINION KLINIPRESS NONWOVEN COMPRESS STERILE 10 X 10 CM 4 PLY PACKED PER 100 X 1 PCS REF 175030 *S*	100 pcs.	1600 pcs.			175030
3017054	KLINION KLINIPRESS NONWOVEN COMPRESS STERILE 10 X 10 CM 4 PLY PACKED PER 50 X 2 PCS REF 175031 *S*	100 pcs.	1600 pcs.			175031
3017055	KLINION KLINIPRESS NONWOVEN COMPRESS STERILE 10 X 10 CM 4 PLY PACKED PER 20 X 5 PCS REF 175033 *S*	100 pcs.	1600 pcs.			175033
3016481	KLINION KLINIPRESS NON WOVEN COMPRESS STERILE 10 X 10 CM PACKED PER 5 PCS 4 PLY REF 40512072 *S*	250 pcs.	2000 pcs.			40512072
3017076	KLINION KLINIPRESS NONWOVEN COMPRESS EXTRA STERILE 10 X 10 CM 8 PLY PACKED PER 100 X 1 PCS REF 175130 *S*	100 pcs.	1600 pcs.			175130
3017077	KLINION KLINIPRESS NONWOVEN COMPRESS EXTRA STERILE 10 X 10 CM 8 PLY PACKED PER 50 X 2 PCS REF 175131 *S*	100 pcs.	1600 pcs.			175131
3017058	KLINION KLINIPRESS NONWOVEN COMPRESS STERILE 10 X 20 CM 4 PLY PACKED PER 50 X 1 PCS REF 175040 *S*	50 pcs.	600 pcs.			175040
3017059	KLINION KLINIPRESS NONWOVEN COMPRESS STERILE 10 X 20 CM 4 PLY PACKED PER 25 X 2 PCS REF 175043 *S*	50 pcs.	600 pcs.			175041
3017079	KLINION KLINIPRESS NONWOVEN COMPRESS EXTRA STERILE 10 X 20 CM 8 PLY PACKED PER 50 X 1 PCS REF 175140 *S*	50 pcs.	600 pcs.			175140
3017062	KLINION KLINIPRESS DRAIN COMPRESS NONWOVEN 5 X 5 CM 4 PLY PACKED PER 50 X 2 PCS REF 175050 *S*	100 pcs.	2000 pcs.			175050
3017063	KLINIPRESS NW DRAINCOMPRES SPLITKOMPRES 10X10CM 4LAGEN STERIEL 50X2ST 175051	100 pcs.	1600 pcs.			175051