

## KLINION - Barrier Cream

### Description


Moisturizing, liquid-repellent cream with 20% fat content for protection of dry and sensitive skin. The Klinion Personal Care Barrier liquid-repellent Cream can be used on dry and sensitive skin for protection and/or prevention of redness. The lotion is particularly well suited for protection of visible areas of the skin exposed to wetting, such as around the mouth and nose. Environmentally friendly and non-allergenic ingredients make Barrier Cream very suitable for use in children and infant care.

### Product information




In compliance with: Cosmetic product regulation EC 1223/2009

#### Features

- Nordic Swan Eco-label
- Perfume-free
- Available in 3 different sizes; 15, 100 and 150 ml
- REF: 4100 is a barrier cream with milk serum



## Article overview

Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
 3015688	KLINION PERSONAL CARE BARRIER CREAM PERFUME FREE 15ML 30808413	30 pcs.	30 pcs.			30808413
 3015687	KLINION PERSONAL CARE BARRIER CREAM PERFUME FREE 100ML 30808412	12 pcs.	12 pcs.			30808412
 3016607	KLINION PERSONAL CARE BARRIER CREAM TUBE MILK SERUM 150 ML REF 4100	1 pcs.	6 pcs.			4100