

KLINION - Wash Cream

Description

Wash cream for cleaning the skin without water. pH 5,5. The Klinion Personal Care Body Wash Cream is a good alternative to soap and water if a patient is bedridden as it is not necessary to rinse the cream off after application. A suitable amount should be placed on either a dry or moist single use cloth or napkin and massaged into the skin. Alternatively, it can be applied directly onto the skin. There is no need to apply further moisturizing agents onto the skin after using Klinion Wash Cream.

Product information


In compliance with: Cosmetic product regulation EC 1223/2009

Features

- Nordic Swan Eco-label
- Perfume-free
- Available in 600 ml



Article overview

Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
 3015684	KLINION PERSONAL CARE WASH CREAM FOR SKIN CLEANSING PH 5.5 600ML 30808409	10 pcs.	10 pcs.			30808409