

KLINIPRESS - Novopad

Description

Nonwoven wound compress. For superficial wounds with little to moderate exudate. With a non-adhering film layer on both sides. Fixate with bandage or an adhesive securing plaster. Medical device. Sterile if pouch and seal are intact.

Product information

In compliance with: MDD (EU) 93/42/EEC Class Is



Article overview

	Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
	3017064	KLINION KLINIPRESS NOVOPAD NON ADHERENT DRESSING 5 X 5 CM PACKED PER 100 X 1 PCS 175060 *S*	100 pcs.	1600 pcs.			175060
	3017065	KLINION KLINIPRESS NOVOPAD NON ADHERENT DRESSING 10 X 10 CM PACKED PER 100 X 1 PCS 175061 *S*	100 pcs.	1600 pcs.			175061
	3017066	KLINION KLINIPRESS NOVOPAD NON ADHERENT DRESSING 10 X 20 CM PACKED PER 50 X 1 PCS REF 175062 *S*	50 pcs.	600 pcs.			175062