

## MESITRAN - Honey products

### Description

L-Mesitran is a pioneering product range, based on Medical Grade Honey, available worldwide since 2002 with excellent clinical results.

These advanced wound care dressings from Holland were the first to receive the CE mark and FDA approval. They have been used on a variety of severe wounds around the globe, including antibiotic resistant (MRSA) infected wounds. L-Mesitran currently is incorporated in many standard wound care protocols around the world.

L-Mesitran offers a natural breakthrough in wound care.

### Product information

Six types of products can be used to treat all wounds. Please find a protocol instruction below.

In compliance with: MDD (EU) 93/42/EEC Class IIb



## Article overview

Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
3017022	KLINION ADVANCED L-MESITRAN WONDZALF 50GR REF 174700	1 pcs.	12 pcs.			174700
3036571	KLINION ADVANCED L-MESITRAN WONDZALF 20GR REF 174701	1 pcs.	14 pcs.			174701
3017026	KLINION ADVANCED L-MESITRAN SOFT HONINGGEL 50GR REF 174706	1 pcs.	12 pcs.			174706
3036549	KLINION ADVANCED L-MESITRAN SOFT HONINGGEL 15GR REF 174705	1 pcs.	14 pcs.			174705
3017027	KLINION ADVANCED L-MESITRAN HYDRO HONINGVERBAND 10X10CM REF 174710 *S*	5 pcs.	100 pcs.			174710
3017028	KLINION ADVANCED L-MESITRAN BORDER HONINGVERBAND MET FIXATIERAND 10X10CM REF 174711 *S*	5 pcs.	100 pcs.			174711
3017029	KLINION ADVANCED L-MESITRAN NET HONINGNETVERBAND 10X10CM REF 174712 *S*	5 pcs.	100 pcs.			174712
3017030	KLINION ADVANCED L-MESITRAN TULLE ZALFKOMPRES 10X10CM REF 174715	10 pcs.	240 pcs.			174715