


KLINIDERM - Superabsorbent Non-Sterile

Description

Absorbent wound dressing for heavily exuding wounds; effective moisture uptake and distribution, comfortable, does not adhere to the wound, does not leak due to special moisture-repellent outer layer. White = wound side; green = external side. Do not cut or tear the dressing. Medical device.

Product information

In compliance with: MDD (EU) 93/42/EEC Class IIb



Article overview

	Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
	3016456	KLINION KLINIDERM SUPERABSORBENT DRESSING NON STERILE 10 X 10 CM REF 40511710	100 pcs.	600 pcs.			40511710
	3016457	KLINION KLINIDERM SUPERABSORBENT DRESSING NON STERILE 10 X 20 CM REF 40511711	50 pcs.	300 pcs.			40511711
	3016458	KLINION KLINIDERM SUPERABSORBENT DRESSING NON STERILE 20 X 20 CM REF 40511712	50 pcs.	150 pcs.			40511712