


## KLINIDERM - Foam

### Description

Foam dressing for the treatment of chronic and acute wounds. The foam dressing is soft, comfortable and highly absorbent. Creates the ideal environment for moist wound healing and avoids adherence to the wound surface. Latex-free. Medical device. Sterile if pouch and seal are intact.

### Product information

In compliance with: MDD (EU) 93/42/EEC Class IIb



## Article overview

Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
3051275	KLINIDERM FOAM 5X5CM 40174810 STERIEL	10 pcs.	100 pcs.			40174810
3051276	KLINIDERM FOAM 10X10CM 40174811 STERIEL	10 pcs.	100 pcs.			40174811
3051277	KLINIDERM FOAM WOUND DRESSING 10X18CM REF 40174813 STERILE	10 pcs.	100 pcs.			40174813
3051278	KLINIDERM FOAM 15X15CM 40174812 STERIEL	10 pcs.	100 pcs.			40174812
3051293	KLINION ADVANCED KLINIDERM FOAM SCHUIMVERBAND HIEL 13X19.5CM 40174820 STERIEL	5 pcs.	40 pcs.			40174820
3051294	KLINIDERM FOAM 20X20 40174814 STERIEL	10 pcs.	100 pcs.			40174814