



Datasheet and IFU

REF 9151313

Content	Packaging	Size
Iris scissor straight, sharp/sharp	1/25/100	11,5



Storage:

Dry and cool without any exposure to UV light.

IFU:

The scissor has a sharp end and a handle. The scissor is used for cutting sutures, swabs etc.

If there is any sign of damages on the packaging the goods must be discarded

Intended use: cutting sutures and other kind of material and only to be used by professionals

Packaging	qty	Material	Disposal
Primary	1	Plastic /medical paper	combustion/hard plastic
Secondary	25	TDF carton	carton
Tertiary	4x25 = 100	TCF carton	carton
Art no		9151313	
Manufacturer		Shaoxing Yibon Medical Ltd, Kina	
Notified Body		0197 TÜV, Rheinland (Shanghai) CO.,Ltd	
CE mark acc		MDD 93/42/ECC annex IX	
Notified Body		0197 TÜV, Rheinland (Shanghai) CO.,Ltd	
Sterilization		EO acc ISO 11135:2014	
Sterile		yes	
Expiry		5 years from production	
Content of latex/pvc		No/no	
Reusable		No, single use	

Symbols used

Medical device	
EO Sterilization acc ISO 11135:2014	
Single use only	
Ce mark acc MDD 93/42/EEC class IS. Notified body 0197 TÜV Rheinland (shanghai) Co., ltd	
Sterilization	
The packaging is partially reusable	

Manufacturer:
Shaoxing Yibon Medical Co., Ltd
Paojiang Industrial Zone
Shaoxing, Zhejiang, China

Importer:
Danomed ApS
Rørvangsvej 51
4300 Holbæk, Denmark

EC Rep
Obelis S.A
BD. Général Wahis 53
1030 Brussels, Belgium