


Mixing Set

- Different lengths of tubing between 180-205 cm
- Roller clamp intended to ensure consistent and safe infusion rate
- Spike without silicone to enable safe connection to infusion bag/bottle/vial
- Liquid filter and spike with air filter
- Closed male luer adapter to enable safe preparation, administration and disconnection **
- Flowstop with hydrophobic filter, designed to simplify priming and reduce the risk of leakage and contamination
- Self-ejecting rotating luer-lock **
- MRI compatible
- 20 drp aqua dest = 1 ml

Mixing Set, non-vented

Art no.	Description	PVC-free	Ward/Trsp.box
60MR03508057	Mixing Set 180 cm, tube to vial 20 cm		10/100*
60MR03508057L	Mixing Set 180 cm, tube to vial 35 cm		10/100*



Art no.	Description	DEHP-free PVC	Ward/Trsp.box
60MR03508357T	Mixing Set 180 cm, tube to vial 25 cm		10/100*

