

sharpsafe®

The all new 7L Sharps Container

With DUAL PROTECT **PRO™**

DATA SHEET



The Sharpsafe® 7L container is ideally suited for sharps disposal in medium use areas and ideal for acute settings. Offering industry leading features such as the new Dual Protect Pro™, the container is designed to encourage horizontal fill to maximise capacity. The new container design has a 13% greater capacity than the previous model. Combined with the low carbon footprint recycled material in the base, this container supports carbon reduction goals within healthcare settings.

Lid colour options:

Yellow

Product Code: 51125430 GTIN Code: 5056271508986

Orange

Product Code: 51125410 GTIN Code: 5056271508993

Purple

Product Code: 51125420 GTIN Code: 5056271509006

Blue

Product Code: 51125440 GTIN Code: 5056271509310

Blue Scottish

Product Code: 51125220 GTIN Code: 5056271509921

Red Scottish

Product Code: 51125280 GTIN Code: 5056271509938

Key Features:

- ✓ The first sharps container to feature Dual Protect Pro™ 1
- ✓ In-Mould Label
- ✓ 100% Recycled Plastic Base
- ✓ Luer Lock/Luer Slip Disconnection Feature 2
- ✓ Star Disconnection Feature 3
- ✓ Integrated Handle
- ✓ Customisable Label (on request)
- ✓ Four Click Assembly

Safer today for a sustainable tomorrow.



Vernacare
clinical waste
management

Pocket Small Medium Large Specialist Accessories

The all new 7L Sharps Container

Certifications

| | |
|------------------|-----|
| ISO 23907-1:2019 | Yes |
| BSI Kitemark | Yes |
| UN3291 | Yes |
| ISO 9001 | Yes |
| ISO 14001 | Yes |
| ISO 45001 | Yes |
| Medical Device | No |
| CE Marked | N/A |

Conformities

| | |
|----------|---------|
| AS 23907 | Pending |
| TRBA 250 | Yes |

Logistical Data

| | |
|----------------------------------|----------|
| Carton Qty | 48 |
| Weight When Empty (kg) | 0.84 |
| Net Weight (kg) | 20.4 |
| Carton Width (mm) | 356 |
| Carton Height (mm) | 712 |
| Carton Length (mm) | 595 |
| Volume (m ³) | 0.151 |
| MOQ | 1 Carton |
| Lead Time | 10 Days |
| Storage Temperature Conditions | Ambient |
| UK Pallet – Cartons Per Pallet | 10 |
| UK Pallet – Units Per Pallet | 480 |
| Euro Pallet – Cartons Per Pallet | 8 |
| Euro Pallet – Units Per Pallet | 384 |

Product Information

| | |
|----------------------------|---------|
| Total Volume (L) | 6.96 |
| Usable Volume (L) | 5.85 |
| Capacity (%) | 84 |
| Weight of Base (g) | 292.4 |
| Weight of Lid (g) | 113 |
| Combined Weight (g) | 405.4 |
| Fill Line Height (mm) | 287 |
| Height (mm) | 341.9 |
| Width (mm) | 170 |
| Depth (mm) | 205 |
| Aperture Size (mm) (L x W) | 61 x 65 |
| Sterile | No |
| Lot/Batched | Yes |
| Expiry Date | N/A |

Material Information

| | |
|---------------------------|------------------------|
| Base Material | Recycled Polypropylene |
| Lid Material | Polypropylene |
| Latex Free | Yes |
| Free From Animal Products | Yes |
| REACH Statement | Yes |

Accessory Compatibility

| Accessory | Product Code |
|------------------|--------------|
| Stand | 42570000 |
| SafeClip Bracket | 45000000 |
| Wall Loop | Coming Soon |

Pocket
Small
Medium
Large
Specialist
Accessories

Safer today for a sustainable tomorrow.

1 Western Avenue, Matrix Park, Buckshaw Village, Chorley, PR7 7NB, UK
+44 (0)1204 529494 hello@vernacare.com vernacare.com

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice. © Copyright Vernacare International Limited 2024. All rights reserved. CWM009 REV03

