


SUCTION SET

REF. NO. **2000 1800 10 ML US**



COMPONENT DESCRIPTION

COMPONENT NO.

1. SUCTION PIECE

Yankauer "DuCatheter" 4 relieving eyes, 2 x bent

0.00025095

OD: 9.3 mm, length: 250 mm

Material: PVC (GC1202D9N2)


2. FUNNEL 9.7

COMPONENT NO.

Material: FEA 508 N T141 – PVC 50°ShoreA, DEHP-FREE

0.00000097

Dye material: Eva No. 161666, blue



REQUIREMENT (mm)

A	B	C	D	E	F	G
16.8 (± 0.1)	13.9 (± 0.1)	13.3 (± 0.1)	12.7 (± 0.1)	9.53 (± 0.1)	8.98 (± 0.1)	7.0 (± 0.1)
H	I	J	K	L	M	
46.0 (± 0.1)	10.0 (± 0.1)	6.0 (± 0.1)	9.7 (± 0.1)	13.8 (± 0.1)	8.7 (± 0.1)	

3. TUBE


COMPONENT NO.

MultiLumen transparent tube.

Material: FEA 768 N T01 – PVC 76° ShoreA, DEHP – FREE,

- 2000 1800 10 ML US – length: 1700 (± 17) mm

0.70170010

	
Ø d INSIDE DIAMETER	Ø D OUTSIDE DIAMETER
7.0 (- 0.2) mm	10.0 (- 0.3) mm

PACKAGING MATERIAL

COMPONENT NO.

Inner bag 165 x 550 mm

2.16555004

Label on the inner bag 35 x 70 mm

7.00003570

Outer bag 310 x 600 mm

25 pcs/bag

2.31060010

Box 386 x 286 x 196 mm

25 pcs/box

6.38628619

Label on the box 85 x 131 mm

7.00085131

Label US on the box 34 x 15 mm

7.00003415

Work instruction

QI 01

Packaging instruction

QI 04

Comments:

Funnel and tube

Glued

Not glued

X

Suction piece and funnel

X

Number of coils:

1800 mm – **2 coils.**

Suction pieces (rigid tube) put separately on tube.

	TITLE: DEVICE MASTER RECORD	
	QUALITY FILE	Ref. no.: QF 04-5
	Revision: 3	Page 3 of 3
	Issue date: 24.03.2025	Effective date: 07.11.2025

Prepared by:	Joanna Kaczmarek	Signature:	
Reviewed by:	Beata Krasowska	Signature:	
Approved by:	Monika Łochyńska	Signature:	