


SUCTION SET

REF. NO. **2000 0900 10 ML US**



COMPONENT DESCRIPTION

COMPONENT NO.

1. SUCTION PIECE

Yankauer "DuCatheter" 4 relieving eyes, 2 x bent

0.00025095

OD: 9.3 mm, length: 250 mm

Material: PVC (GC1202D9N2)


2. FUNNEL 9.7

COMPONENT NO.

Material: FEA 508 N T141 – PVC 50°ShoreA, DEHP-FREE

0.00000097

Dye material: Eva No. 161666, blue



	TITLE: DEVICE MASTER RECORD	
	QUALITY FILE	Ref. no.: QF 04-5
	Revision: 1	Page 2 of 2
	Issue date: 22.04.2025	Effective date: 22.04.2025

3. TUBE


COMPONENT NO.

MultiLumen transparent tube.

Material: FEA 768 N T01 – PVC 76° ShoreA, DEHP – FREE,

- 2000 0900 10 ML US – length: 850 (± 8.5) mm

0.70085010

	
$\varnothing d$ INSIDE DIAMETER	$\varnothing D$ OUTSIDE DIAMETER
7.0 (- 0.2) mm	10.0 (- 0.3) mm

PACKAGING MATERIAL

COMPONENT NO.

Inner bag 150 x 360 mm

2.15036004

Label on the inner bag 35 x 70 mm

7.00003570

Outer bag 310 x 600 mm

25 pcs/bag

2.31060010

Box 386 x 286 x 196 mm

25 pcs/box

6.38628619

Label on the box 85 x 131 mm

7.00085131

Label US on the box 34 x 15 mm

7.00003415

Work instruction

QI 01

Packaging instruction

QI 04

Comments:

	Glued	Not glued
Funnel and tube	X	
Suction piece and funnel		X
Suction pieces (rigid tube) put separately on tube.		

Prepared by:	Joanna Kaczmarek	Signature:	
Reviewed by:	Beata Krasowska	Signature:	
Approved by:	Monika Łochyńska	Signature:	