

# Product Specifications

Brand : Klinion®  
 Product : Kliniplast® Plaster Strips

### Product description

Water-resistant wound plaster strips. Contains four different plaster sizes.

### Intended use

For protecting wounds and wound surfaces.



### Classification

Class I under the Medical Devices Regulations (SI 2002 No. 618, as amended) and  
 Class I under Regulation (EU) 2017/745 Medical Device Regulation (MDR).

### Product & Packaging Reference

Product Reference	Description	Size (*)	Piece(s)	Inner Carton (pcs)	Outer Carton (pcs)
40294114	KLINIPLAST® PLASTER STRIP ASSORTI	72 x 19 MM	4	20	720
		72 x 25 MM	8		
		38 x 9 MM	6		
		Ø 22MM	2		

\*Sizes are approximate

Product Reference	Description	Size (*)	EAN Inner Carton	EAN Outer Carton
40294114	KLINIPLAST® PLASTER STRIP ASSORTI	72 x 19 MM	8717438370040	8717438370057
		72 x 25 MM		
		38 x 9 MM		
		Ø 22 MM		

\*Sizes are approximate

# Product Specifications

## Features

- Water-resistant
- Breathable
- Elastic
- Adheres well
- Easy to remove
- No significant adhesive residue

## Precautions

- Single use
- Carefully clean wound and surrounding skin.
- Dry the skin area before applying the dressing
- Replace when necessary, at least daily

## Specifications

Component	Material
Backing material	Plastic PE backing
Adhesive	Acrylic
Wound pad	Nonwoven pad covered with non-adherent film: 70% PP. 30% Viscose
Release liner	Paper liner with pull tab
Wrapper	Individually sealed paper wrapper

Performance Characteristic	72 x 19 MM	72 x 25 MM	38 x 9 MM	Ø 22 MM	Specification	Test method
Tensile strength	>0.7 kg/19mm	>0.9kg/25mm	>0.3kg/9.5mm	N/A	Kg/mm	ISO 29862
Adhesive strength	600-1500 g/19mm	800-2000 g/25mm	300-750 g/9.5mm	N/A	g/mm	EN 13726-3

# Product Specifications

## Presence of specific components

Component	Present
Latex	No
Phthalates	No
Substances of animal origin	No
Substances of human origin	No
Blood or blood derivatives	No
Medicinal substances	No
Nanomaterials	No
Heavy metals	No
Colorants/Dyestuffs	No
Silicone	No

## Sterilization

N.A

## Storage conditions

The products should be stored in the original packaging, in a cool, well-ventilated and dry place, protected from direct sunlight.

## Shelf life

Kliniplast® Plaster Strips has a shelf-life of 5 years from production date, when stored under specified conditions.

Approval	Name	Date	Signature
Product manager wound care	Florin Lazar	<b>Florin Lazar</b>	 Digitally signed by Florin Lazar Date: 2025.09.30 15:13:13 +02'00'
RA specialist	Emil Laske	<b>Emil Laske</b>	 Digitally signed by Emil Laske Date: 2025.09.30 14:08:53 +02'00'