



Altius High Pressure Central Venous Catheter



ALTIUS HP					Priming Volume (ml)					Flow Rates (ml/hr)				
Product Code	Description	Lumens	French size	Length (cm)	Proximal	Distal	Medial 1	Medial 2	Medial 3	Proximal	Distal	Medial 1	Medial 2	Medial 3
					K2CS115501P	ALTIUS HP 1 Lumen CVC Set 5Fr 15cm	1	5	15		0.3			
K2CS120501P	ALTIUS HP 1 Lumen CVC Set 5Fr 20cm	1	5	20		0.3					2800(16G)			
K2CS115701P	ALTIUS HP 1 Lumen CVC Set 7Fr 15cm	1	7	15		0.55					9300(14G)			
K2CS120701P	ALTIUS HP 1 Lumen CVC Set 7Fr 20cm	1	7	20		0.6					7500(14G)			
K2CS115702P	ALTIUS HP 2 Lumen CVC Set 7Fr 15cm	2	7	15	0.3	0.35				2500(17G)	5800(14G)			
K2CS120702P	ALTIUS HP 2 Lumen CVC Set 7Fr 20cm	2	7	20	0.3	0.35				2050(17G)	5100(14G)			
K2CS115802P	ALTIUS HP 2 Lumen CVC Set 8Fr 15cm	2	8	15	0.3	0.4				3200(16G)	6600(14G)			
K2CS120802P	ALTIUS HP 2 Lumen CVC Set 8Fr 20cm	2	8	20	0.3	0.4				2700(16G)	6100(14G)			
K2CS115703P	ALTIUS HP 3 Lumen CVC Set 7Fr 15cm	3	7	15	0.25	0.3	0.25			1600(18G)	3200(16G)	1600(18G)		
K2CS120703P	ALTIUS HP 3 Lumen CVC Set 7Fr 20cm	3	7	20	0.25	0.3	0.25			1350(18G)	2800(16G)	1350(18G)		
K2CS130703P	ALTIUS HP 3 Lumen CVC Set 7Fr 30cm	3	7	30	0.25	0.3	0.25			1100(18G)	2350(16G)	1100(18G)		
K2CS115803P	ALTIUS HP 3 Lumen CVC Set 8Fr 15cm	3	8	15	0.25	0.25	0.3			1600(18G)	6600(14G)	1600(18G)		
K2CS120803P	ALTIUS HP 3 Lumen CVC Set 8Fr 20cm	3	8	20	0.25	0.25	0.3			1350(18G)	5100(14G)	1350(18G)		
K2CS115904P	ALTIUS HP 4 Lumen CVC Set 9Fr 15cm	4	9	15	0.25	0.25	0.25	0.25		1600(17G)	6600(14G)	1600(17G)	1600(17G)	
K2CS120904P	ALTIUS HP 4 Lumen CVC Set 9Fr 20cm	4	9	20	0.25	0.25	0.25	0.25		1350(17G)	6050(14G)	1350(17G)	1350(17G)	
K2CS115905P	ALTIUS HP 5 Lumen CVC Set 9Fr 15cm	5	9	15	0.2	0.3	0.2	0.2	0.2	1250(18G)	3200(16G)	1250(18G)	1250(18G)	1250(18G)
K2CS120905P	ALTIUS HP 5 Lumen CVC Set 9Fr 20cm	5	9	20	0.2	0.3	0.2	0.2	0.2	1100(18G)	2800(16G)	1100(18G)	1100(18G)	1100(18G)
K2CS130905P	ALTIUS HP 5 Lumen CVC Set 9Fr 30cm	5	9	30	0.2	0.3	0.2	0.2	0.2	840(18G)	2200(16G)	840(18G)	840(18G)	840(18G)



Altius High Flow Central Venous Catheter



ALTIUS HP HF

Product Code	Description	Lumens	French size	Length (cm)	Priming Volume (ml)					Flow Rates (ml/hr)				
					Proximal	Distal	Medial 1	Medial 2	Medial 3	Proximal	Distal	Medial 1	Medial 2	Medial 3
K2CS1151255PHF	ALTIUS HP HF 5 Lumen CVC Set 12.5Fr 15cm	5	12.5	15	0.35	0.38	0.48	0.48	0.38	1600(18G)	3200(16G)	11400(12G)	11400(12G)	3200(16G)
K2CS1201255PHF	ALTIUS HP HF 5 Lumen CVC Set 12.5Fr 20cm	5	12.5	20	0.42	0.45	0.5	0.5	0.45	1350(18G)	2750(16G)	10500(12G)	10500(12G)	2750(16G)
K2CS1251255PHF	ALTIUS HP HF 5 Lumen CVC Set 12.5Fr 25cm	5	12.5	25	0.5	0.48	0.55	0.55	0.48	1100(18G)	2610(16G)	10200(12G)	10200(12G)	2610(16G)

ALTIUS HF

Product Code	Description	Lumens	French size	Length (cm)	Priming Volume (ml)					Flow Rates (ml/hr)				
					Proximal	Distal	Medial 1	Medial 2	Medial 3	Proximal	Distal	Medial 1	Medial 2	Medial 3
K2CS115853HF	ALTIUS HF 3 Lumen CVC Set 8.5Fr 15cm	3	8.5	15	0.38	0.48	0.38			3200(16G)	3200(16G)	11400(12G)		
K2CS120853HF	ALTIUS HF 3 Lumen CVC Set 8.5Fr 20cm	3	8.5	20	0.45	0.5	0.45			2750(16G)	2750(16G)	10500(12G)		
K2CS125853HF	ALTIUS HF 3 Lumen CVC Set 8.5Fr 25cm	3	8.5	25	0.48	0.55	0.48			2610(16G)	2610(16G)	10200(12G)		

For more information, pricing or samples, speak to your local Kimal Representative, or contact:

© 2025 Kimal PLC. All rights reserved. No part of this document may be reproduced in any form or by any means without permission in writing from Kimal PLC.
Kimal PLC reserve the right to change the specification of its products at any time.
Kimal PLC, Clayfield Road, Worcester Six Business Park, Worcestershire, WR4 0AE.

UK
Tel: +44 (0)845 437 9542
Fax: +44 (0)845 437 9544
Email: customercare@kimal.com

International
Tel: +44 (0)845 437 9543
Fax: +44 (0)845 437 9549
Email: internationalenquiries@kimal.com