

DELTA VEN® NEO DELTA VEN®

Peripheral Intravenous Catheter of FEP and PUR



 **DELTA MED**
MEDICAL DEVICES
Solutions for Health


PART OF



NUMANTEC
HUMAN INNOVATION TECHNOLOGY

PERIPHERAL IV CATHETER

DELTA VEN® NEODELTA VEN®



The I.V. Catheters division manufactures high quality IV cannulae and medical devices guaranteeing best safety and comfort for both patients and healthcare professionals.


Available in 3 configurations



Straight




Straight & winged




Ported & winged

Catheter material:

- clear
 - flexible
 - anti kinking
 - visible under X-Ray (3 radiopaque lines)
- Available in Polyurethane (PUR) version (particularly heat-sensitive: softening once in vein) and Teflon (FEP) version.



Back-cut, sharp needle bevel for smooth, non-traumatic insertion



Translucent hub for easy flashback visualization

Main Benefits:

- Delta Ven and Neo Delta Ven are Power Injectable (no 26G)
- Wide range of gauge-lengths and configurations

PERIPHERAL IV CATHETER



Neo Delta Ven 2 / Delta Ven 2 - Ported and winged intravenous catheter

NEO DELTA VEN 2	DELTA VEN 2	DESCRIPTION	GAUGE	CATHETER LENGHT (mm) (inch)	CATHETER OD (mm)	CATHETER ID (mm)	FLOW RATE (ml/min)	SALES UNIT / UNIT PER BOX
3104022	3103222	Ported and winged	26	19 mm (3/4")	0,64	0,45	19	50/200
3113122	3111522	Ported and winged	24	19 mm (3/4")	0,74	0,55	29	50/200
3123122	3121522	Ported and winged	22	25 mm (1")	0,9	0,65	42	50/200
3133122	3131522	Ported and winged	20	32 mm (1" 1/4)	1	0,75	59	50/200
3183122	3181522	Ported and winged	18	32 mm (1" 1/4)	1,3	0,95	103	50/200
3143122	3141522	Ported and winged	18	45 mm (1" 3/4)	1,3	0,95	96	50/200
3153122	3151522	Ported and winged	17	45 mm (1" 3/4)	1,5	1,15	155	50/200
3163122	3161522	Ported and winged	16	45 mm (1" 3/4)	1,75	1,35	225	50/200
3173122	3171522	Ported and winged	14	45 mm (1" 3/4)	2	1,55	290	50/200



Neo Delta Ven 1 / Delta Ven 1 - Straight winged catheter

NEO DELTA VEN 1	DELTA VEN 1	DESCRIPTION	GAUGE	CATHETER LENGHT (mm) (inch)	CATHETER OD (mm)	CATHETER ID (mm)	FLOW RATE (ml/min)	SALES UNIT / UNIT PER BOX
3108122	3106522	Straight and winged	26	19 mm (3/4")	0,64	0,45	19	50/200
3118122	3116522	Straight and winged	24	19 mm (3/4")	0,74	0,55	29	50/200
3128122	3126522	Straight and winged	22	25 mm (1")	0,9	0,65	42	50/200
3138122	3136522	Straight and winged	20	32 mm (1" 1/4)	1	0,75	59	50/200
3188122	3186522	Straight and winged	18	32 mm (1" 1/4)	1,3	0,95	103	50/200
3148122	3146522	Straight and winged	18	45 mm (1" 3/4)	1,3	0,95	96	50/200
3168122	3166522	Straight and winged	16	45 mm (1" 3/4)	1,75	1,35	225	50/200
3178122	3176522	Straight and winged	14	45 mm (1" 3/4)	2	1,55	290	50/200



Neo Delta Ven T / Delta Ven T - Straight catheter

NEO DELTA VEN T	DELTA VEN T	DESCRIPTION	GAUGE	CATHETER LENGHT (mm) (inch)	CATHETER OD (mm)	CATHETER ID (mm)	FLOW RATE (ml/min)	SALES UNIT / UNIT PER BOX
3108522	3106422	Straight	26	19 mm (3/4")	0,64	0,45	19	50/200
3118522	3116422	Straight	24	19 mm (3/4")	0,74	0,55	29	50/200
3128522	3126422	Straight	22	25 mm (1")	0,9	0,65	42	50/200
3138522	3136422	Straight	20	32 mm (1" 1/4)	1	0,75	59	50/200
3188522	3186422	Straight	18	32 mm (1" 1/4)	1,3	0,95	103	50/200
3148522	3146422	Straight	18	45 mm (1" 3/4)	1,3	0,95	96	50/200
3168522	3166422	Straight	16	45 mm (1" 3/4)	1,75	1,35	225	50/200
3178522	3176422	Straight	14	45 mm (1" 3/4)	2	1,55	290	50/200

MKIVC00362 - 17/09/2024 - REV. 4.0

Delta Med
is a **Numantec**
company

DELTA MED
MEDICAL DEVICES
Solutions for Health

medical supplies and accessories
Securmed

PentaFerte
MEDICAL DEVICES SOLUTIONS

D.B.M.
INNOVATIVE
MEDICAL DEVICES

adriamed
EXCELLENCE IN MEDICAL DEVICES

PHS MEDICAL



NUMANTEC[®]
HUMAN INNOVATION TECHNOLOGY

CERTIFICATES:
ISO 13485 - ISO 9001 - MDSAP 13485
The EC Certificates for Medical Devices manufactured by
Delta Med are issued by TÜV SÜD Product Service GmbH.

Delta Med SpA
via Guido Rossa, 20 46019 Viadana (MN), Italy
Telefono: +39 0375 785915 - Fax: +39 0375 785201
info.deltamed@numantec.eu - www.deltamed.it
info@numantec.eu - www.numantec.eu