


PRODUKTDATABLAD

SUG MEDENA M0300 SUGESLANGE 2MX5,5X8MM M/NIPPEL



SUG MEDENA M0300 SUGESLANGE 2MX5,5X8MM M/NIPPEL

Sug Medena M0300 sugeslange 2mx5,5x8mm med nippel. Medena sugeslange 2 meters lengde, diameter 5,5x8mm. Nippel uten fingerkontrollhull.

Varenr.	Tekst	Enhet
702516	Sug Medena M0300 sugeslange 2mx5,5x8mm m/nippel	100

MEDIQ NORGE AS

Dyrskuevegen 13, Pb. 113, N-2041 Kløfta
Tlf: (+47) 67 02 43 00 - E-mail: kundeservice.no@mediq.com
www.mediqnorge.no