

## SPECIFICATION

### Pouch Rolls for Hot-Air Sterilisation 100 m x 100 mm



---

This pouch cannot be sealed using a sealing machine, so it is recommended that it be sealed using the appropriate sealing tape prior to hot-air sterilisation. Flat type, without gusset, for temperatures up to 180° C.

#### IMPORTANT NOTE:

Only for use with hot-air sterilisers in the cosmetic field. IT IS NOT A MEDICAL DEVICE!

## SPECIFICATION

### Pouch Rolls for Hot-Air Sterilisation

Dimensions	PZN	REF	Unit
100 m x 30 mm		F5 36030	1 roll
100 m x 50 mm		F5 36050	1 roll
100 m x 75 mm		F5 36075	1 roll
100 m x 100 mm		F5 36100	1 roll
100 m x 150 mm		F5 36150	1 roll
100 m x 200 mm		F5 36200	1 roll
100 m x 250 mm		F5 36250	1 roll

