

KLINIPRESS - HG Gauze Compress Sterile

Description

Gauze swab for cleansing and covering. Highly absorbent. 100% cotton. Change daily. Medical device. Sterile if pouch and seal are intact.















Product information



In compliance with: MDD (EU) 93/42/EEC Class Is



Article overview

	Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
	3016682	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 5 X 5 CM 8 PLY PACKED PER 50 X 2 PCS REF 111010 *S*	100 pcs.	3000 pcs.			111010
	3016487	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 5 X 5 CM PACKED PER 2 PCS 8 PLY REF 40512090 *S*	160 pcs.	4800 pcs.			40512090
	3016684	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 5 X 5 CM 8 PLY PACKED PER 30 X 5 PCS REF 111015 *S*	150 pcs.	3000 pcs.			111015
	3016685	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 5 X 5 CM 12 PLY PACKED PER 10 X 1 PCS REF 111016 *S*	10 pcs.	200 pcs.			111016
	3016686	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 5 X 5 CM 12 PLY PACKED PER 100 X 1 PCS REF 111017 *S*	100 pcs.	2000 pcs.			111017
	3016687	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 5 X 5 CM 12 PLY PACKED PER 50 X 2 PCS REF 111018 *S*	100 pcs.	2000 pcs.			111018
	3016689	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 5 X 5 CM 12 PLY PACKED PER 20 X 5 PCS REF 111020 *S*	100 pcs.	1200 pcs.			111020
	3016699	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 7.50 X 7.50 CM 8 PLY PACKED PER 50 X 2 PCS REF 111041 *S*	100 pcs.	2000 pcs.			111041
	3016488	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 7.50 X 7.50 CM PACKED PER 2 PCS 8 PLY REF 40512091 *S*	200 pcs.	3000 pcs.			40512091

	Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
	3016701	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 7.50 X 7.50 CM 8 PLY PACKED PER 30 X 5 PCS REF 111043 *S*	150 pcs.	2250 pcs.			111043
	3016492	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 7.50 X 7.50 CM PACKED PER 5 PCS 8 PLY REF 40512101 *S*	200 pcs.	3000 pcs.			40512101
	3016691	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 7.50 X 7.50 CM 12 PLY PACKED PER 20 X 5 PCS REF 111022 *S*	100 pcs.	4000 pcs.			111022
	3016679	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 8.50 X 5 CM 8 PLY PACKED PER 16 X 1 PCS REF 111001 *S*	16 pcs.	384 pcs.			111001
	3016680	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 8.50 X 5 CM 8 PLY PACKED PER 100 X 1 PCS REF 111002 *S*	100 pcs.	1600 pcs.			111002
	3016705	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 10 X 10 CM 8 PLY PACKED PER 100 X 1 PCS REF 111082 *S*	100 pcs.	1600 pcs.			111082
	3016704	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 10 X 10 CM 8 PLY PACKED PER 50 X 2 PCS REF 111080 *S*	100 pcs.	2000 pcs.			111080
	3016489	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 10 X 10 CM PACKED PER 2 PCS 8 PLY REF 40512092 *S*	200 pcs.	1600 pcs.			40512092
	3016702	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 10 X 10 CM 8 PLY PACKED PER 30 X 5 PCS REF 111078 *S*	150 pcs.	1500 pcs.			111078
	3016706	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 10 X 10 CM 12 PLY PACKED PER 25 X 1 PCS REF 111084 *S*	25 pcs.	400 pcs.			111084
	3016707	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 10 X 10 CM 12 PLY PACKED PER 100 X 1 PCS REF 111085 *S*	100 pcs.	1600 pcs.			111085
	3016708	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 10 X 10 CM 12 PLY PACKED PER 50 X 2 PCS REF 111086 *S*	100 pcs.	1600 pcs.			111086
	3016710	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 10 X 10 CM 12 PLY PACKED PER 20 X 5 PCS REF 111088 *S*	100 pcs.	1200 pcs.			111088
	3016711	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 10 X 10 CM 12 PLY PACKED PER 10 X 10 PCS REF 111089 *S*	100 pcs.	800 pcs.			111089
	3016715	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 10 X 20 CM 12 PLY PACKED PER 50 X 1 PCS REF 111112 *S*	50 pcs.	600 pcs.			111112
	3016716	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 10 X 20 CM 12 PLY PACKED PER 40 X 2 PCS REF 111113 *S*	80 pcs.	1280 pcs.			111113

	Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
	3016725	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 20 X 20 CM 8 PLY PACKED PER 10 X 5 PCS REF 111235 *S*	50 pcs.	800 pcs.			111235
	3016724	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 20 X 20 CM 16 PLY PACKED PER 30 X 1 PCS REF 111223 *S*	30 pcs.	480 pcs.			111223
	3016726	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 20 X 25 CM 16 PLY PACKED PER 25 X 1 PCS REF 111255 *S*	25 pcs.	400 pcs.			111255
	3016727	KLINION KLINIPRESS HG GAUZE COMPRESS STERILE 20 X 30 CM 16 PLY PACKED PER 20 X 1 PCS REF 111256 *S*	20 pcs.	160 pcs.			111256
	3016982	KLINION KLINIPRESS GAUZE DRAIN COMPRESS 10 X 10 CM 16 PLY PACKED PER 1 PCS REF 170115 *S*	50 pcs.	600 pcs.			170115
	3016983	KLINION KLINIPRESS GAUZE DRAIN COMPRESS 10 X 10 CM 16 PLY PACKED PER 2 PCS REF 170116 *S*	50 pcs.	600 pcs.			170116