

## KLINION - Sharps Containers

### Description

The Klinion Easy Care Sharps containers are designed for hospitals and health care institutions to dispose of syringes, scalpels and other types of sharps.

### Product information

#### Features

- Design and manufacturing secure a high degree of safety
- Products are made of strong and durable material
- Containers are waterproof when closed
- Containers can be sterilized



### Article overview

	Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
	3028327	KLINION EASY CARE SHARPS CONTAINER 0.20 LITER YELLOW 77510900	1 pcs.	100 pcs.			77510900
	3028328	KLINION EASY CARE SHARPS CONTAINER 0.50 LITER YELLOW 77510905	1 pcs.	100 pcs.			77510905
	3028329	KLINION EASY CARE SHARPS CONTAINER 1.50 LITER YELLOW 77510910	1 pcs.	100 pcs.			77510910
	3028330	KLINION EASY CARE SHARPS CONTAINER 2 LITER YELLOW 77510915	1 pcs.	100 pcs.			77510915
	3028332	KLINION EASY CARE SHARPS CONTAINER 5 LITER YELLOW 77510921	1 pcs.	30 pcs.			77510921
	3028334	KLINION EASY CARE SHARPS CONTAINER 7 LITER YELLOW 77510926	1 pcs.	25 pcs.			77510926
	3028336	KLINION EASY CARE SHARPS CONTAINER 11 LITER YELLOW 77510931	1 pcs.	20 pcs.			77510931
	3028338	KLINION EASY CARE SHARPS CONTAINER 21 LITER YELLOW 77510936	10 pcs.	10 pcs.			77510936