

## KLINIDERM - Film Dressings

### Description

Transparent wound film for surface wounds with slight exsudate. Creates an ideal wound climate. Transparent, semi-permeable and self-adhesive. Fixation (I.V.). Medical device. Sterile if pouch and seal are intact.

### Product information

Item 3017013 will be replaced with 3052512: Wound film for fixation of intravenous cannula (IV). Semi-permeable, self-adhesive, with transparent window and nonwoven edges. Sterile if pouch and seal are intact.

In compliance with: MDD (EU) 93/42/EEC Class Is



## Article overview

Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
3016547	KLINION KLINIDERM FILM DRESSING TRANSPARENT STERILE 4 X 5 CM REF 40514840 *S*	50 pcs.	600 pcs.			40514840
3016548	KLINION KLINIDERM FILM DRESSING TRANSPARENT STERILE 5 X 7.20 CM REF 40514841 *S*	50 pcs.	600 pcs.			40514841
3016549	KLINION KLINIDERM FILM DRESSING TRANSPARENT STERILE 10 X 10 CM REF 40514842 *S*	50 pcs.	600 pcs.			40514842
3016550	KLINION KLINIDERM FILM DRESSING TRANSPARENT STERILE 10 X 12 CM REF 40514843 *S*	50 pcs.	600 pcs.			40514843
3016551	KLINION KLINIDERM FILM DRESSING TRANSPARENT STERILE 10 X 15 CM REF 40514844 *S*	50 pcs.	600 pcs.			40514844
3016552	KLINION KLINIDERM FILM DRESSING TRANSPARENT STERILE 10 X 20 CM REF 40514845 *S*	50 pcs.	300 pcs.			40514845
 3016553	KLINION KLINIDERM FILM DRESSING TRANSPARENT STERILE 10 X 25 CM REF 40514846 *S*	50 pcs.	300 pcs.			40514846
 3016554	KLINION KLINIDERM FILM DRESSING TRANSPARENT STERILE 10 X 30 CM REF 40514847 *S*	50 pcs.	300 pcs.			40514847
3016555	KLINION KLINIDERM FILM DRESSING TRANSPARENT STERILE 15 X 20 CM REF 40514848 *S*	50 pcs.	300 pcs.			40514848

Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
3016556	KLINION KLINIDERM FILM DRESSING TRANSPARENT STERILE 20 X 30 CM REF 40514849 *S*	50 pcs.	300 pcs.			40514849
3016557	KLINION KLINIDERM FILM IV TRANSPARENT FILM INFUSION PLASTER 6 X 7 CM REF 40514855 *S*	50 pcs.	600 pcs.			40514855
3052512	KLINION KLINIDERM FILM IV NON WOVEN INFUUSPLEISTER TRANSPARANTE FOLIE MET VERSTEVIG 7.5X8.9CM 174420 STERIEL	50 pcs.	300 pcs.			40174420