

## KLINIDUR - Short Stretch

### Description

Short stretch bandage. Ambulatory compression therapy, provides high pressure during activity and low pressure at rest. Medical device.

### Product information

In compliance with: Regulation (EU) 2017/745 Medical Device Regulation (MDR) Class I



## Article overview

Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
3058770	KLINION KLINIDUR TEXTELAST COMPRESSION BANDAGE SHORT STRETCH SKIN COLOR 5 M X 8 CM PACKED PER 2 PCS REF 132628	2 pcs.	50 pcs.			132628
3058771	KLINION KLINIDUR TEXTELAST COMPRESSION BANDAGE SHORT STRETCH COLOR BEIGE 5 M X 10 CM PACKED PER 2 PCS REF 132630	2 pcs.	42 pcs.			132630
3016596	KLINION KLINIDUR COMPRESSION BANDAGE SHORT STRETCH 5 M X 8 CM REF 40516002	10 pcs.	120 pcs.			40516002
3016597	KLINION KLINIDUR COMPRESSION BANDAGE SHORT STRETCH 5 M X 10 CM REF 40516003	10 pcs.	80 pcs.			40516003
3016598	KLINION KLINIDUR COMPRESSION BANDAGE SHORT STRETCH 5 M X 12 CM REF 40516004	10 pcs.	80 pcs.			40516004
3067593	KLINION KLINIDUR COHESIVE COMPRESSION BANDAGE 6CMX5M REF 40511802 NON- STERILE	1 pcs.	24 pcs.			40511802
3067594	KLINION KLINIDUR COHESIVE COMPRESSION BANDAGE 8CMX5M REF 40511803 NON- STERILE	1 pcs.	24 pcs.			40511803
3067595	KLINION KLINIDUR COHESIVE COMPRESSION BANDAGE 10CMX5M REF 40511804 NON- STERILE	1 pcs.	24 pcs.			40511804
3067526	KLINION KLINIDUR COHESIVE COMPRESSION BANDAGE 12CMX5M REF 40511805 NON- STERILE	1 pcs.	24 pcs.			40511805