

KLINIGRIP - Arm Sling NonWoven

Description

Triangular bandage, nonwoven, with safety pin. For immobilising arms or shoulders. Medical device.

Product information

In compliance with: Regulation (EU) 2017/745 Medical Device Regulation (MDR) Class I



Article overview

Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
3038155	KLINION KLINIGRIP ARM SLING NON WOVEN (MITELLA) 96 X 96 X 136 CM REF 102026	48 pcs.	288 pcs.			102026