

precisa[®] N shock-proof

Shock-resistant aneroid sphygmomanometer



shock-proof



Handheld, portable, aneroid Sphygmomanometer with shock resistant technology, meeting standard DIN EN 81060-1:2012-08. Validated by the British and Irish Hypertension Society.

- › Shock-resistant according to DIN EN 81060-1:2012-08.
- › Shock-resistant from a height of 1 m (39.37 in).
- › Precision air release valve with fine adjustment, wear-free.
- › Microfilter protects air release valve and measuring system.
- › Diaphragm pressure-loading capacity up to 600 mmHg.
- › Specially tempered copper-beryllium diaphragm with precision movement and almost non-aging.
- › Made in Germany
- › 3-year calibration warranty
- › Latex-free.
- › Metal valve base.
- › Ergonomically shaped spoon-type handle made of stainless steel.
- › Linear tilted aluminium scale Ø 63 mm (2.48 in) from 0 – 300 mmHg, easy to read.
- › Nylon zipper bag

No. 1364-100 precisa[®] N shock-proof w/o cuff

No. 1364-122 precisa[®] N shock-proof, 1-tube, velcro cuff (large adult)

No. 1364-107 precisa[®] N shock-proof, 1-tube, velcro cuff (adult)

No. 1364-192 precisa[®] N shock-proof, 1-tube, velcro cuff (children)

Please visit riester.de for more information