

SPECIFICATION

Mediware Disposable overboot



Protects shoes, socks and trouser legs.

Made of strong, transparent polyethylene. With ties.

Ideal for veterinary medicine and in all areas with strict hygiene requirements.

Shaft height:43 cm

Dispenser box:50 pieces

Carton:500 pieces

SPECIFICATION

Mediware Disposable overboot

Size	PZN	REF	Unit
Universal size		H9 06	50 pcs.

Product	REF	Quantity	Unit price	Total

_____ plus tax

Personal details:

Customer number:

Company:

First name: Last name:

E-mail:

Street, street#:

Postal code, city:

I hereby place a binding order:

Company stamp / Address

Date / Sign

You are ordering under our terms and conditions, our valid terms and conditions can be found at <https://www.servoprax.de/en/terms-and-conditions/>