

SPECIFICATION

DCT rectal catheter



Disposable, with funnel, 2 side eyelets and closed tip.

Product advantages:

- Colour coded
- Ribbed grip
- No lumen reduction at the adapter
- Large, rounded atraumatic eyelet
- Heat-sensitive PVC which conforms to body contours
- Rounded, smooth tip

Length = 40 cm

Sterile and singly packed.

SPECIFICATION**DCT rectal catheter**

Size	OC	PZN	REF	Unit
CH 22	360 pieces	13229230	F5 150622	30 pcs.
CH 25	360 pieces	13229247	F5 150625	30 pcs.
CH 28	360 pieces	13229253	F5 150628	30 pcs.
CH 30	360 pieces	13229282	F5 150630	30 pcs.
CH 32	360 pieces	13229299	F5 150632	30 pcs.

