

412

White colour
Cor branca

Conforming bandages

Ligaduras elásticas

Bandes extensibles

Elastische Fixierbinden

Elastisch fixatiewindsel

Benda elàstica

Venda extensible

Elastisk fixeringsbinda



Mais imagens



Light weight conforming bandages

- Viscose and polyamide.
- For fixation of dressings.
- Construction when stretched: 17 threads/sq.cm.
- Stretched Length: 4m.
- Packing:
 - Non sterile: with or without individual wrapping in cellophane + shelf box + transport carton.
 - Sterile: peel-open pack + shelf box + transport carton.

Non sterile

Code	Size	Pcs./Cellophane	Pcs./Sh.box	Pcs./Tr. carton
412-001	4m x 4cm	1	20	560
412-002	4m x 5cm	1	20	440
412-003	4m x 6cm	1	20	400
412-004	4m x 7cm	1	20	320
412-005	4m x 8cm	1	20	280
412-006	4m x 10cm	1	20	240
412-007	4m x 12cm	1	20	200
412-008	4m x 15cm	1	20	160
412-009	4m x 20cm	1	20	120
412-010	4m x 4cm		20	560
412-011	4m x 5cm		20	440
412-012	4m x 6cm		20	400
412-013	4m x 7cm		20	320
412-014	4m x 8cm		20	280
412-015	4m x 10cm		20	240
412-016	4m x 12cm		20	200
412-017	4m x 15cm		20	160
412-018	4m x 20cm		20	120

Sterile

Code	Size	Pcs./Pack	Pcs./Sh.box	Pcs./Tr. carton
412-501	4m x 4cm	1	40	720
412-502	4m x 5cm	1	40	720
412-503	4m x 6cm	1	45	810
412-504	4m x 7cm	1	150	900
412-505	4m x 8cm	1	150	900
412-506	4m x 10cm	1	130	780
412-507	4m x 12cm	1	110	660
412-508	4m x 15cm	1	100	600
412-509	4m x 20cm	1	80	480