

## KLINION - Impervious Isolation Gowns

### Description

Non-woven, PE-coated impervious isolation gown with knitted cuffs and velcro neck closure, colour yellow. Dispose according to local procedures and guidelines. Not made with natural rubber latex. Medical device.

### Product information

In compliance with: Regulation (EU) 2017/745 Medical Device Regulation (MDR) Class I

Tested in accordance with EN 13795-1:2019



## Article overview

| Article number | Description   | Inner carton | Outer carton | EAN (IC) | EAN (OC) | Reference |
|----------------|---|--------------|--------------|----------|----------|-----------|
| 3027534        | KLINION PROTECTION KLINIPRO IMPERVIOUS GOWN YELLOW M 114 CM REF 741244  | 10 pcs.      | 50 pcs.      |          |          | 741244    |
| 3027533        | KLINION PROTECTION KLINIPRO IMPERVIOUS GOWN YELLOW L 134 CM REF 741243  | 10 pcs.      | 50 pcs.      |          |          | 741243    |
| 3026974        | KLINION PROTECTION KLINIPRO IMPERVIOUS GOWN YELLOW XL 144 CM WITH CUFF AND VELCRO REF 702100                          | 10 pcs.      | 50 pcs.      |          |          | 702100    |
| 3025547        | KLINION PROTECTION ISOLATION APRON LONG SLEEVES WITH TRICOT MANCHET M 115 CM 50 PCS YELLOW WITH PE COATING REF 522240 | 1 pcs.       | 50 pcs.      |          |          | 522240    |
| 3025550        | KLINION PROTECTION ISOLATION APRON LONG SLEEVES WITH TRICOT MANCHET L 135   | 1 pcs.       | 50 pcs.      |          |          | 522250    |