

## KLINION - Aprons

### Description

Disposable apron made of polyethylene (PE) material. Barrier protection against harmless fluids and dirt. Protection against general and minimal risks only, which do not expose user to temperatures > 50°C or to dangerous impacts. Before usage inspect the apron for any defect or imperfections. For single use only. Re-use can cause contamination.











### Product information





In compliance with: PPE Regulation (EU) 2016/425 Category I

- Available in many sizes and thickness
- Also available in short and long sleeve



## Article overview

	Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
	3025536	KLINION PROTECTION APRON PE 110 X 80 CM 100 PCS WHITE REF 522110	100 pcs.	1000 pcs.			522110
	3025537	KLINION PROTECTION APRON PE 125 X 80 CM 100 PCS WHITE REF 522130	100 pcs.	1000 pcs.			522130
	3025538	KLINION PROTECTION APRON PE WHITE 140CM BOX PE 140 X 80 CM 80 PCS WHITE REF 522140	80 pcs.	800 pcs.			522140
	3025539	KLINION PROTECTION APRON PE 165 X 80 WHITE REF 522160	70 pcs.	700 pcs.			522160
	3025565	KLINION PROTECTION PE APRON 71 X 117 CM WHITE EMBOSSED 0.023 MM REF 522700-1	100 pcs.	500 pcs.			522700-1
	3025566	KLINION PROTECTION APRON 81 X 125 CM TRANSPARANT EMBOSSED 0.025 MM REF 522701-1	50 pcs.	500 pcs.			522701-1
	3025739	KLINION PROTECTION PE APRON 117 X 71 CM WHITE 18 MM EMBOSSED 0.018 MM REF 522700	100 pcs.	500 pcs.			522700
	3025740	KLINION PROTECTION APRON 125 X 81 CM TRANSPARANT EMBOSSED 0.0195 MM REF 522701	50 pcs.	500 pcs.			522701
	3025741	KLINION PROTECTION APRON 125 X 81 CM NATURAL EMBOSSED 0.0195 MM REF 522702	50 pcs.	500 pcs.			522702
	3025742	KLINION PROTECTION APRON 160 X 90 CM WHITE EMBOSSED 0.019 MM REF 522703	70 pcs.	420 pcs.			522703

	Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
	3025743	KLINION PROTECTION APRON 140 X 80 CM WHITE FOLDED 0.031 MM REF 522704	70 pcs.	420 pcs.			522704
	3025744	KLINION PROTECTION APRON ON ROLL 150 X 80 CM TRANSPARANT 0.0195 MM REF 522707	50 pcs.	500 pcs.			522707
	3025745	KLINION PROTECTION APRON FOLDED 150 X 80 CM TRANSPARANT 0.0195 MM REF 522708	40 pcs.	400 pcs.			522708
	3049840	KLINION PERSONAL PROTECTION DISPOSABLE LDPE APRON, WHITE 522709	70 pcs.	420 pcs.			522709