

Product specifications

Mediq Norge

Made by Ab Fatex Oy

Quality requirements			206	205
Testing parameter	Dim.	Coc	Sheet TP Regular 80x210cm	Sheet TP Regular 80x140
Appearance, property				
- Color			Blue	Blue
- Dirt			Free from dirt, lint and stains	Free from dirt, lint and stains
- Number of polyester threads			-	-
- Fluid barrier			Both sides 2100mm	Both sides 1400mm
Components				
- Tissue (ECF bleached)			Towel tissue 30 gsm/ 780mm	Towel tissue 30 gsm/ 780mm
- Tissue 2 (ECF bleached)			-	-
- PE-film			PE-film 775 x 0,014mm Blue	PE-film 775 x 0,014mm Blue
- Adhesive (Hotmelt)			2,44g/pcs	1,63g/pcs
- Fluid barrier			0,73g/pcs	0,49g/pcs
- Polyester threads			-	-
Length	mm	Yes	2100mm ± 30mm	1400mm ± 30mm
Width	mm	Yes	780mm ± 10mm	780mm ± 10mm
Length (folded)	mm		525mm	350mm
Width (folded)	mm		195mm	195mm
Weight	g	Yes	66,7g ± 15%	44,5g ± 15%
Water penetration		Yes	0	0
Packaging				
Pcs/Stack (surrounded by banderoll)	Qty		25pcs	25pcs
Stacks/Carton	Qty		4	4
Cartons/Pallet	Qty		64	96
Carton size (Inner)	mm		550x390x130	390x370x130
Banderoll width	mm		47	47
Pathogenic germs		Yes	Free from pathogenic germs like staph. Aureus, Enterobacteriae, E. coli, Pseudomonas aeruginosa	Free from pathogenic germs like staph. Aureus, Enterobacteriae, E. coli, Pseudomonas aeruginosa
Check of certificate		Yes	By request	By request