

## KLINIGRIP - Eco-ideal Bandage

### Description

Universal support bandage. Light support bandage, white. 60% cotton, 39% polyamide, 1% polyurethane, maintains elasticity even after washing (60°C). 5m when stretched. For support following treatment of strains, sprains and injuries. For use when applying a compression bandage in combination with synthetic wadding. Medical device.

### Product information

In compliance with: MDD (EU) 93/42/EEC Class I



## Article overview

Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
3016853	KLINION KLINIGRIP ECO-IDEAL SUPPORT BANDAGE 60% COTTON CREAM COLOR 5 M X 6 CM REF 132469	10 pcs.	40 pcs.			132469
3016854	KLINION KLINIGRIP ECO-IDEAL SUPPORT BANDAGE 60% COTTON WHITE 5 M X 8 CM REF 132470	10 pcs.	40 pcs.			132470
3016855	KLINION KLINIGRIP ECO-IDEAL SUPPORT BANDAGE 60% COTTON WHITE 5 M X 10 CM REF 132471	10 pcs.	40 pcs.			132471
3016856	KLINION KLINIGRIP ECO-IDEAL SUPPORT BANDAGE 60% COTTON WHITE 5 M X 12 CM REF 132472	10 pcs.	40 pcs.			132472