

## KLINIFIX - Cohesive Fixation Bandage

### Description



Cohesive elastic securing bandage, white. Latex free.

### Product information

In compliance with: Regulation (EU) 2017/745 Medical Device Regulation (MDR) Class I



## Article overview

	Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
	3016838	KLINION KLINIFIX COHESIVE ELASTIC FIXATION BANDAGE SELF-ADHESIVE LATEX FREE 4 M X 2,5 CM REF 132377	10 pcs.	450 pcs.			132377
	3016830	KLINIFIX COHESIVE ELASTIC BANDAGE LATEX FREE 4MX4CM 132378	5 pcs.	280 pcs.			132378
	3038996	KLINION KLINIFIX COHESIVE ELASTIC FIXATION BANDAGE LATEX-FREE 4MX6CM REF 132379	5 pcs.	280 pcs.			132379
	3038997	KLINION KLINIFIX COHESIVE ELASTIC FIXATION BANDAGE LATEX FREE 4 M X 8 CM REF 132380	5 pcs.	140 pcs.			132380
	3039000	KLINION KLINIFIX COHESIVE ELASTIC FIXATION BANDAGE LATEX FREE 4 M X 10 CM REF 132381	5 pcs.	140 pcs.			132381
	3039001	KLINIFIX COHESIVE ELASTIC BANDAGE LATEX FREE 4MX12CM 132382	5 pcs.	140 pcs.			132382
	3016835	KLINIFIX COHESIVE ELASTIC BANDAGE LATEX FREE 20MX6CM 132386	1 pcs.	24 pcs.			132386
	3016836	KLINION KLINIFIX COHESIVE ELASTIC FIXATION BANDAGE SELF-ADHESIVE LATEX FREE 20 M X 8 CM REF 132387	1 pcs.	24 pcs.			132387
	3016837	KLINION KLINIFIX COHESIVE ELASTIC FIXATION BANDAGE SELF-ADHESIVE LATEX FREE 20 M X 10 CM REF 132388	1 pcs.	24 pcs.			132388