

SPECIFICATION

Mediware Long stitch cutter long form



Made of stainless steel, long form

Single, sterile in aluminium blister pack

SPECIFICATION

Mediware Long stitch cutter long form

Size	OC	PZN	REF	Unit
Long form, 110 mm	6000 pcs.		I1 0503	100 pcs.

Product	REF	Quantity	Unit price	Total

_____ plus tax

Personal details:

Customer number: _____

Company: _____

First name: Last name: _____

E-mail: _____

Street, street#: _____

Postal code, city: _____

I hereby place a binding order:

Company stamp / Address

Date / Sign

You are ordering under our terms and conditions, our valid terms and conditions can be found at <https://www.servoprax.de/en/terms-and-conditions/>