



Product Data Sheet

for

Thermometer Cover Sheaths

(Article number 403625)

Prepared by:

Minitube AB
March 2023

Product Data Sheet for Thermometer Covers

Product name:	Tempasept with ointment
Product no:	403625
Legal manufacturer:	Minitube AB, Sweden
Fits thermometer types:	Digital
Storing Conditions:	Ambient room temperature (16-25 °C), 36 months

Inner Cover

Cover Material:	Polyethylene film LDPE
Cover Colour:	Transparent
Flaps Material:	Polyethylene film
Flaps Colours:	White, green
Lubricant/ointment:	Yes

Outer Cover

Bottom Material:	Bleached craft paper; clear, water based barrier coating
Bottom Colour:	Blue
Back Film:	Polyester film
Front Material:	Polypropylene film
Front Colour:	Transparent
Front film Imprint:	Yes

Inner Box

Material:	Cardboard
Dimensions:	210 mm x 80 mm x 155 mm
Imprint:	Yes
Contents:	1,000 covers
Weight (empty):	50 g
Weight (filled):	970 g

Outer Box (Skipper)

Material:	Cardboard
Dimensions:	410 mm x 215 mm x 310 mm
Label:	Yes
Weight (empty):	335 g
Weight (filled):	9,7 kg
Contents:	10 inner boxes

Pallets

Quantities:	500,000 Covers
	500 Inner boxes
	50 Outer boxes
Dimensions:	
including EUR pallet:	120 cm x 83 cm x 180cm
excluding pallet:	111 cm x 83 cm x 170 cm
Weight:	
including pallet:	510 kg
excluding pallet:	490 kg