



Sol-M™ Catheter Tip Syringe Without Needle, 3-Piece

Sol-M™ 3-Piece Syringes accurately measure the quantity of liquid medications (according to ISO 7886-1) and are available in a wide range of sizes to address a wide variety of clinical needs.

Features:

- High-transparency barrel to easily visualize contents
- Silicone lubrication ensures smooth plunger movement
- Plunger backstop that prevents accidental withdrawal of plunger from the

barrel

- Large finger grips allow for easy, controlled injections

Intended use: The Catheter Tip Syringe is used to irrigate, withdraw fluid from, or instill fluid into, a body cavity or wound.

Intended user: Licensed healthcare professionals (HCP).

Cautions:

- Single use device. Re-use or use if the package is damaged may lead to infection or other illness/ injury.

CE 2797

[Product Listing](#)

[Videos](#)

[Catalog](#)

[IFU](#)

Ref. Number	Description	Box Qty	Case Qty
180060CT	Sol-M Catheter Tip Syringe	30	240
18100CT	Sol-M 100 ml catheter Tip Syringe	30	60