




TITLE: DEVICE MASTER RECORD	
QUALITY FILE	Ref. no.: QF 04-5
Revision: 1	Page 1 of 1
Issue date: 13.07.2020	Effective date: 13.07.2020

ORTHOPEDIC SUCTION SET MED EXTRA FILTER

REF. NO. 1832 3000 80 ML / 1832 3500 80 ML



COMPONENT DESCRIPTION

COMPONENT NO

1. ORTHOPEDIC SUCTION KIT

- a. 2 filters- material: SAN
- b. 1 handle- material: SAN
- c. 3 tips- PVC :

i. ID/OD: 6.0/8.0 mm (CH24) 150 mm	0.06015080
ii. ID/OD: 6.0/8.0 mm (CH24) 230 mm	0.06023080
iii. ID/OD: 6.0/8.0 mm (CH24) 280 mm	0.06028080

2. TUBE

Flexible transparent tube with ribs ID/OD 5.6/8.0 mm, material: PVC 73°ShoreA , DEHP – FREE,

- 1832 3000 80 ML – 2900 mm length 0.56290080
- 1832 3500 80 ML – 3400 mm length 0.56340080

3. DISTAL END/PROXIMAL END 2x funnel 8.0 PVC

material: PCV 50°ShoreA, DEHP – FREE , blue 0.00000080

PACKAGING MATERIALS

Innerbag 250 x 500 mm	2.25050004
Peelpack 270 x 520 mm	5.27052412
Label on the peelpack 85 x 131 mm	7.00085131
Outer bag 310 x 700 mm	2.31070010
Box 590 x 386 x 260 mm	6.59038626
Label on the box 85 x 131 mm	7.00085131

12 or 13 pcs/bag
25 pcs/box

Operating instructions	QI 01
Packing instructions	QI 04