

NO. 20195-034

Bipolar forceps Classic, bayonet

Bipolar forceps Classic, bayonet, tip 2 mm, blunt, length 190 mm



Technical characteristics

Packaging unit	1 pcs.
Cleaning in washer disinfector	max 93 °C
Sterilization in autoclave	max 134 °C
Shape / Forming	Bayonet
Observe instructions for use / manual	AB
Max. electrical capacity	0,5 kVp
Total length of instrument	190 MM
Tip	Tip 2 mm, blunt
Gripping length	6 MM