

Skater™ Fix


Reduce risk of catheter dislocation – without suturing

SKATER™ FIX is a latex free catheter fixation device, consisting of a central thin carrier which has hydrocolloid adhesive and is covered with non-woven medical grade adhesive. The soft foam support block and velcro prevent the catheter from kinking.

The catheter will be mechanically fixed with velcro, and minimal tape is involved directly on the catheter.

The velcro cover can be opened for puncture site inspection. In this way you have the possibility to look for infection or leakage whenever required.

- Latex free
- Safe anchoring of drainage catheter up to one week
- Allows inspection of puncture site without change of dressing
- Easy to change
- Minimal risk of catheter dislocation during change of fixation disk
- Foam bridge prevents catheter kinking
- Clean catheter with no tape residuals
- Adjustable to catheters from 5 to 16 French
- Hydrocolloid carrier will absorb body fluids from puncture site
- Minimum risk of skin irritation
- Disc follows patient contour and movements



DRAINAGE

SKATER™ FIX catheter fixation device
(Box of 25)

CATALOGUE NO.	DESCRIPTION
802500001	Latex free catheter fixation device

Skin Fix™ for Connecting Tubes

Skin Fix™ consists of an adhesive area of microporous tape which is applied to the patient's skin, and a fixation strap which keeps the connecting tubes in place.

- Comfortable (non-irritant)
- Permeable to humidity and heat
- Moisture resistant
- Follows the contours of the skin
- Fixation strap – with strong adhesion



Skin Fix™ (Box of 25)

CATALOGUE NO.	DESCRIPTION
802500000	Sterile, single packed