



Product description

Order number	85.1019
Product description	Safety lancet neonatal, pink, blade 1.5 mm, penetration depth 1.2 mm - disposable latest generation penetration aid for collecting capillary blood
PZN number	06796164

Product characteristics

Example of use	Children, Adults with large order quantities, neonates: small blood count, blood sugar, bilirubin
Puncture site	Heel, finger
Penetration depth	1,2 mm

Size

Blood volume	50-250 µl
Width of product	21 mm
Length of product	52 mm
Width of blade	1,5 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product in accordance with MDD or IVDD	Medical Device
CE certified	CE0124
Sterilisation	sterile
Lot or serial number	batched

Packaging

Minimum order qty.	1200
Pcs. / inner box	200
Pcs. / box	1200
Pcs. / pallet	108000
Depth of case	366 mm
Width of case	248 mm
Height of case	173 mm
Case volume	0,0157 cbm
Weight of case	4,08 KG
EAN of inner box	4038917085069
EAN of outer box	4038917057660

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

