

Drainage

Connecting Tubes, Adapters, and Stopcocks

We offer a wide range of connecting tubes for all types of drainage bags. The connecting tubes are available with or without stopcocks. Flexible, kink-resistant tubing; clinical versatility.

Box of 25

Catalogue #	Description	Catalogue #	Description
401012010	10cm Female-Male Luer Lock	404040015	15cm male - male luer-lock
401012025	25cm Female -Male Luer Lock	404040030	30cm male - male luer-lock
401012050	50cm Female - Male Luer Lock	404043015	15cm male-3-way stopcock/male luer-lock
401012080	80cm Female - Male Luer Lock	404043030	30cm male-3-way stopcock/male luer-lock
401013025	25cm 1-way stopcock/Female-Male	404343015	15cm male luer-lock /3-way stopcock - funnel
401013050	50cm 1-way stopcock/ Female-Male	404343030	30cm male luer-lock /3-way stopcock - funnel
401013080	80cm 1-way stopcock/ Female-Male	404743015	15cm male luer-lock - olive
403943015	15cm female-olive	80700103	3-way stopcock - low pressure, max pressure 50 psi
404017015	15cm male-funnel		
404017030	30cm male-funnel		
404018015	15cm one-way stopcock/male-funnel		
404018030	30cm one-way stopcock/male-funnel		
404022000	Male/three-way stopcock/funnel 6cm		
404024000	Male/funnel 6cm		

Connecting Tubes - Low Pressure (Box of 10)

Catalogue #	Description
802020000	Femaleolive with stopper
802021000	Maleolive with stopper

