


800430 OP-Towel

Product details

Size : 45cm x 75cm
Descriptive feature : 2Ply, Absorbent
Sterile : Sterile
Brand : BARRIER

Images



Delivered items

800430-60

Country of origin: Thailand

Shelf life: 5 years

Sterilization method: Beta

Packing information: First packaging layer is a peel-open sterile barrier, paper/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a corrugated board dispenser box. Third layer is a corrugated board transport box.

Is suitable for Tray: No

Packing level	Quantity	GS1 code	Width x Length x Height	Vol	Weight gross / net
Piece	1	7323190257782			
Consumer pack	1	7323190272594			
Dispenser box	25	7323190262113			
Transport box	150	7323190262106	298x683x203 mm	41.3 dm3	4.9 / 3.8 kg

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2021 Mölnlycke Health Care AB. All rights reserved.



Packing level	Quantity	GS1 code	Width x Length x Height	Vol	Weight gross / net
Pallet	4800	7323190262090	800x1200x1768 mm	1697.3 dm3	183.1 / 120.0 kg

Material

Animal tissues :	No
Human blood derivatives :	No
Natural rubber latex :	No
Medicinal substances :	No
Phthalates :	No
Polyvinyl chloride :	No

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2021 Mölnlycke Health Care AB. All rights reserved.

