



Product description

| | |
|---------------------|--|
| Order number | 55.476.005 |
| Product description | Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PS, round base, transparent, push cap, 1,000 piece(s)/bag |

Product characteristics

| | |
|--------------------------|------------|
| Graduation | no |
| Closure type | push cap |
| Base shape | round base |
| Volume of work | 5 ml |
| Centrifugation max (RCF) | 4000 x g |

Size

| | |
|-------------------|-------|
| Diameter | 12 mm |
| Length of product | 75 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|------------------|
| Product material | Polystyrene (PS) |
| Colour of product | transparent |

Purity & certification

| | |
|---------------------------|--------------------------|
| Satisfies the requirement | sealed at 0.6 bar/10 min |
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 2000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 1000 |
| Piece (s) / outer case | 2000 |
| Piece(s) / pallet | 84000 |
| Depth of case | 398 mm |
| Width of case | 398 mm |
| Height of case | 250 mm |
| Case volume | 0.0396 cbm |
| Weight of product | 0.0028 kg |
| Weight of case | 6.19 kg |
| EAN of inner box | 4038917107310 |
| EAN case | 4038917011778 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

