



Product description

Order number	77.580
Product description	Urine container 2.0 l, including loose label enclosed

Product characteristics

Label/ Print	With paper label
Colour of print/label	Yellow/black
Graduation	Yes
Type of cap	Screw cap
Light protection	Yes
Design	German/English label
Application	Urine collection
Example of use	Collection and storage of urine
Viewing window	No

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Sample volume	2000 ml
Max. volume	2 l
Diameter of opening	70 mm
Width of product	106 mm
Height of product	240 mm
Length of product	106 mm
Length including cap	238 mm
Length excluding cap	235 mm

Material & colours

Product material	Polyethylene (PE)
Colour of product	Brown
Cap material	Polypropylene (PP)
Colour of cap	gelb

Purity & certification

Product in accordance with MDD or IVDD	IVD
CE certified	CE
Lot or serial number	non batched

Packaging

Minimum order qty.	30
Pcs. / inner box	30
Pcs. / box	30
Pcs. / pallet	360
Depth of case	563 mm
Width of case	343 mm
Height of case	500 mm
Case volume	0,0966 cbm
Weight of case	3,5 KG
EAN of inner box	4038917082655
EAN of outer box	4038917001700

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

