

GLORIOUS ENTERPRISES

PRODUCT SHEET

Products Name: Dressing Forceps, 14.5 cm, serrated



Glorious Art nr. 44.0340.14

Finish: Satin

Net weight: 0.032 gram (approximate)

Material: Stainless steel AISI 410

Sterility: Non-Sterile (Re-Usable)

Exposure Temperature: 134°C (274°F)

Exposure Time: 3 Minutes

CE Classification: Class-I

Storage: 16°C to 30°C (Store the instruments in a closed cabinet, protected from dust, moisture and temperature fluctuations)

life-span: 20 times of use:- after each use, product need to sterilize for ready to use in next operation or next surgery.

Disposal: No special decomposition or disposal required for Surgical Instruments as they don't contain any toxic or hazardous material