


# Product Data Sheet

Bakterie- och virusfilter engångs

REF 28350



Tryckt: 7.4.2014

REF OneMed 287873																	
Egenskaper:	<p>Vitalographs bakterie och virusfilter förhindrar att spirometer eller annan lungfunktionsutrustning kontamineras av patientens utandning, vilket medför mindre risk för smittspridning.</p> <h2>Bacterial Viral Filters</h2> <h3>Technical Data</h3> <table border="1"><tr><td>Product</td><td>BV Filter</td></tr><tr><td>Model Number</td><td>2820</td></tr><tr><td>Efficiency</td><td>99.99% at 55 and 750L/min</td></tr><tr><td>Validation using Bacillus Subtilis var niger mean dia.</td><td>0.7 to 1 µm</td></tr><tr><td>Resistance to Airflow at 750L/min (BVF alone)</td><td>0.08kPa/L/s (0.8 cmH<sub>2</sub>O)</td></tr><tr><td>Resistance to Airflow at 750L/min (BVF + Fleisch):</td><td>0.12kPa/L/s (1.2 cmH<sub>2</sub>O)</td></tr><tr><td>Effective Deadspace (excluding Tapers)</td><td>65ml</td></tr><tr><td>Material</td><td>Electrostatic nanofibre technology filter medium</td></tr></table>	Product	BV Filter	Model Number	2820	Efficiency	99.99% at 55 and 750L/min	Validation using Bacillus Subtilis var niger mean dia.	0.7 to 1 µm	Resistance to Airflow at 750L/min (BVF alone)	0.08kPa/L/s (0.8 cmH <sub>2</sub> O)	Resistance to Airflow at 750L/min (BVF + Fleisch):	0.12kPa/L/s (1.2 cmH <sub>2</sub> O)	Effective Deadspace (excluding Tapers)	65ml	Material	Electrostatic nanofibre technology filter medium
Product	BV Filter																
Model Number	2820																
Efficiency	99.99% at 55 and 750L/min																
Validation using Bacillus Subtilis var niger mean dia.	0.7 to 1 µm																
Resistance to Airflow at 750L/min (BVF alone)	0.08kPa/L/s (0.8 cmH <sub>2</sub> O)																
Resistance to Airflow at 750L/min (BVF + Fleisch):	0.12kPa/L/s (1.2 cmH <sub>2</sub> O)																
Effective Deadspace (excluding Tapers)	65ml																
Material	Electrostatic nanofibre technology filter medium																
Steril produkt:																	
Hållbarhetstid:																	
MBA/st:	50																
Avdelningsförpackning/st:	50																
Transportlåda/st:																	