

GLORIOUS ENTERPRISES

PRODUCT SHEET

Products Name: PEAN Artery Forceps 24cm straight, Smooth jaws GQ



Glorious Art nr. 43.0622.24GQ
Net weight: 110 gram (approximate)
Material: Stainless steel AISI 410
Finish: SAND BLAST
Sterility: Non-Sterile (Re-Usable)
Exposure Temperature: 134°C (274°F)

Exposure Time: 3 Minutes

CE Classification: Class-I

Storage: 16°C to 30°C (Store the instruments in a closed cabinet, protected from dust, moisture and temperature fluctuations)

Shelf life: 3 years

life-span: 20 times of use:- after each use, product need to sterilize for ready to use in next operation or next surgery.

Disposal: No special decomposition or disposal required for Surgical Instruments as they don't contain any toxic or hazardous material