

GLORIOUS ENTERPRISES

PRODUCT SHEET

Products Name: Mosquito Forceps 12.5cm STRAIGHT GQ



Glorious Art nr.	43.0494.12
Finish:	Satin
Net weight:	0.025 gram (approximate)
Material:	Stainless steel AISI 410
Sterility:	Non-Sterile (Re-Usable)
Exposure Temperature:	134°C (274°F)
Exposure Time:	3 Minutes
CE Classification:	Class-I
Storage:	16°C to 30°C (Store the instruments in a closed cabinet, protected from dust, moisture and temperature fluctuations)
life-span:	20 times of use:- after each use, product need to sterilize for ready to use in next operation or next surgery.
Disposal:	No special decomposition or disposal required for Surgical Instruments as they don't contain any toxic or hazardous material