



# **GLORIOUS ENTERPRISES**

## **PRODUCT SHEET**

**Products Name:** COTTON DRESSING JAR WITH LOCK 10X10CM



<b>Glorious Art nr.</b>	<b>60.819.00</b>
<b>Net weight:</b>	<b>205 gram (approximate)</b>
<b>Material:</b>	<b>Stainless steel AISI 202</b>
<b>Finish:</b>	<b>BRIGHT</b>
<b>Sterility:</b>	<b>Non-Sterile (Re-Usable)</b>
<b>Exposure Temperature:</b>	<b>134°C (274°F)</b>
<b>Exposure Time:</b>	<b>3 Minutes</b>
<b>CE Classification:</b>	<b>Class-I</b>
<b>Storage:</b>	<b>16°C to 30°C</b> (Store the instruments in a closed cabinet, protected from dust, moisture and temperature fluctuations)
<b>Shelf life:</b>	<b>3 years</b>
<b>life-span:</b>	<b>20 times of use:- after each use, product need to sterilize for ready to use in next operation or next surgery.</b>
<b>Disposal:</b>	<b>No special decomposition or disposal required for Surgical Instruments as they don't contain any toxic or hazardous material</b>



**P.O. Box 63, Naverland 11, DK 2600 Glostrup, Denmark**  
**Tel: +45 43 62 26 93, Fax: +45 43 62 26 11, Email: [glorious@glorious-enterprises.dk](mailto:glorious@glorious-enterprises.dk)**