



Excellence through research

mondomed[®]


HIGH TECH MEDICAL PVA FOAM

Mondomed RETRACTOR[®] Pad

RETRACTOR
Cover



ABDOMINAL SURGERY



SUPPORTING
MEDICAL
PROFESSIONALS
IN TAKING CARE
OF PEOPLE...
since 1983

PIONEERING KNOW-HOW IN PVA

- MONDOMED® is specialized in the production of high quality polyvinyl alcohol (PVA) foam products for medical applications and has more than 30 years of experience.
- MONDOMED's continuous research and attention to quality has led to excellence and reliability. Today, MONDOMED® is considered to be the only PVA producer in the world offering such a wide range of high tech PVA foam products for medical applications.



Excellence through research
mondomed®
HIGH TECH MEDICAL PVA FOAM

MONDOMED® PVA

BENEFITS

MONDOMED® polyvinyl alcohol (PVA) has unique technical and physical properties which makes the material very suitable for a wide range of medical applications, such as foam devices for abdominal surgery.

MONDOMED® PVA foam products:

- have a spongy, porous structure with a small pore size
- minimize tissue ingrowth
- have an impressive fluid absorption and retaining capacity
- have a smooth foam texture and are atraumatic

MONDOMED® PVA foam is:

- dust and fiber-free
- biocompatible
- environmentally friendly (biodegradable)



Mondomed RETRACTORPad®

MONDOMED® PVA ABDOMINAL SURGERY CUSHION


INDICATIONS

The Retractorpad is a MONDOMED® PVA cushion for abdominal surgery designed for optimal exposure of the surgical area whilst providing excellent fluid absorption and retention and optimal protection of the intestines and blood vessels during abdominal and pelvic surgery. Intestines are kept out of the field of surgery and as such working conditions are facilitated for the surgical team.

PROPERTIES

The Retractorpad is constructed from polyvinyl alcohol in the form of an absorbent microporous foam. As a result its shape is adapted to fit the internal configuration of an abdominal transverse section. Its U-shaped structure allows the Retractorpad to support on the paracolic gutters while the arched device prevents the underlying blood vessels from being compressed. Because of the consistency of the material, there is a good support. As such the intestines will not appear in the field of surgery during the operation but they will stay localised outside this field of surgery. The Retractorpad is slightly moist and flexible which prevents damages to the intestines. Due to its spongy structure, it can absorb blood and other bodily fluids. Also dehydration is avoided and so are contact traumas. If necessary, the Retractorpad can be cut to size, as the material is not fibrous.

- Excellent fluid absorption and retention
- Shaped for optimal intestine support
- U-shaped design prevents compression on the underlying blood vessels
- Soft and moist, to avoid dehydration and contact trauma
- MONDOMED® PVA (Polyvinyl Alcohol) foam is lint and fiber free and can easily be cut to-size
- Offers a great overview of the surgical area, thus facilitating the task of the operating team
- Environmentally friendly



Mondomed RETRACTORPad®

MONDOMED® PVA ABDOMINAL SURGERY CUSHION


Great FLUID ABSORPTION

High FLUID RETENTION capacity

EASY FIT to the abdominal contour (pliable); can be cut to size (fiber free)

OPTIMAL FORM: no pressure on the retro-peritoneal vessels or indirect pressure on the diaphragm

VERY PRACTICAL: optimal protection of the intestines, optimal working conditions for the surgeons



INDICATIONS

- A sterile, sock-like sleeve made out of PVA intended to cover the point and the blade of a surgical retractor during surgery in order to prevent damages to the tissues and organs (to reduce tissue trauma/tissue adhesion during retraction of tissues)

PROPERTIES

- The RetractorCover need to be moistened prior to use.
- After moistening it becomes flexible and has a soft texture which prevents damages to the intestines.
- Because of its porous structure, it can absorb blood and other bodily fluids.



Retractorpad

- 1355910 S
(1 piece/pouch, 25 pouches/box)
- 1355911 L
(1 piece/pouch, 15 pouches/box)



The Retractorpad is

- Pre-moisturised
- Sterilized by irradiation
- X-Ray detectable due to barium thread
- Double packed in two peel pouches

RetractorCover

- 1155200 25x180 mm
- 1155205 38x180 mm
- 1155210 50x180 mm
- 1155211 50 x 240 mm
- 1155215 80x180 mm
- 1155220 65x235 mm
- 1155225 35x185 mm (rounded)
- 1155231 Intestinal Clamp Cover (6 x 90 mm)
- 1155232 18 x 80 mm
- 1155291 58x280 mm

all RetractorCovers are packed:

- 1 piece/pouch
- 10 pouches/folding box
- 9 folding boxes/transport box

Intestinal Clamp Covers are packed:

- 1 pair/pouch
- 10 pouches/folding box
- 9 folding boxes/transport box

The RetractorCover is

- Sterilized by irradiation
- X-Ray detectable due to barium thread
- Double packed in two peel pouches






Mondomed RETRACTORPad[®]

RETRACTOR_{Cover}

Official Distributor



Excellence through research

mondomed[®]

HIGH TECH MEDICAL PVA FOAM



mondomed nv
Middenweg 20
3930 HAMONT-ACHEL
BELGIUM
0032 11 44 09 90
info@mondomed.be
www.mondomed.be

The informations, icons, photographs, images, texts and other in this catalog are protected by industrial and/or intellectual rights and are the property of MONDOMED[®]. As such, any reproduction, representation, adaptation, translation and / or modification in whole or partial is prohibited.
2020 All rights reserved MONDOMED[®]



2012853 08/04/2020

Produced in Belgium

