

microlife®

A partner for people. For life.



Irregular heartbeat detection

IHB

Early alerts to possible heart irregularity.



Product video



BP B2 Basic Blood Pressure Monitor

- IHB technology for detection of irregular heartbeat
- Optimal speed and pressure control with Gentle+
- Comfortable cuff for most arm sizes (22–42 cm)
- 30 data memory and average-display
- Cuff position check

Key features



Irregular heartbeat detection for early alerts to possible heart irregularity



Gentle+ ensures optimal speed and pressure control for comfortable measurements



Easy to read classification of blood pressure result



Washable cuff cover



30 data memory and average-display of all stored readings



Cuff fit check



Cuff check indicator



Gentle+

- Optimal speed and pressure control
- Prevents overpumping for highest comfort during blood pressure measurement



Soft Cuff: 22–42 cm

- Conical shape for better fit and optimal accuracy
- Washable cuff cover for longlasting, hygienic use



Clinically Validated

Highest measuring accuracy with clinical A/A grade according to the British and Irish Hypertension Society (BIHS).

12+
years



Validated and recommended for children. Also validated and recommended for patients with diabetes, renal disease, pregnancy / pre-eclampsia and obesity.



5 year
guarantee