

## Cuffs and supplies

Cuffs, as important as the blood pressure monitor



### Rigid cuff

Conical, preshaped upper arm cuff for high level of comfort

#### Features

- Wide-range cuff size
- Conical, fan-shaped cuff for optimal fit
- Compatible with all Microlife blood pressure monitors
- Rigid material for easy self-attachment
- Clinically validated
- Available size:  
M-L (22–42 cm/8.7–16.5")

[www.microlife.com/rigid-cuff](http://www.microlife.com/rigid-cuff)

### Soft cuff

Conical, soft upper arm cuff with washable cuff cover

#### Features

- Washable cuff cover for long-term and hygienic use
- Conical, fan-shaped cuff for optimal fit
- Compatible with all Microlife blood pressure monitors
- Metal bar for easy self-attachment
- Clinically validated in all sizes
- Different sizes available:  
S (17–22 cm/6.7–8.7")  
M (22–32 cm/8.7–12.6")  
M-L (22–42 cm/8.7–16.5")  
L (32–42 cm/12.6–16.5")  
L-XL (32–52 cm/12.6–20.5")

[www.microlife.com/soft-cuff](http://www.microlife.com/soft-cuff)